



# MICROSCOOP<sup>®</sup> / MINT

The First Subcellular "Pickable" Microscope



## TECHNICAL SPECIFICATIONS

DESCRIPTION		SPECIFICATIONS
<b>Function</b>		Optoproteomics: Ultra-content high-speed microscopy-guided subcellular photoaffinity labeling for hypothesis-free high-precision proteomic discovery
<b>Workflow</b>		Cyclic procedure of the following step: 1. Microscopy imaging: image acquisition 2. Pattern segmentation: selection of user-defined regions of interest 3. Patterned scanning illumination: point-by-point photochemical reactions 4. Stage movement: change of the field of view
<b>Components</b>		<ol style="list-style-type: none"> <li>1 Microscop<sup>®</sup> system (optical engine and electrical controller)</li> <li>2 Inverted epifluorescence microscope Filter sets for microscope</li> <li>3 Epifluorescence illumination light source</li> <li>4 Two-photon laser for Microscop<sup>®</sup> photolabeling</li> <li>5 Camera</li> <li>6 Software package</li> </ol>
1	<b>Dimensions (L x W x H)</b>	Electrical controller : 44 cm x 22 cm x 47 cm Optical engine: 68 cm x 46 cm x 22 cm
	<b>Power Source</b>	100 - 240 VAC, 50/60 Hz
2	<b>Objectives</b>	10x (up to NA 0.45) 20x (up to NA 0.80) 40x (up to NA 0.95)
	<b>Stage</b>	Motorized XY positioning stage (X: ±57 mm, Y: ±36.5 mm stroke) with a vessel holder, suitable for microscope slides, chamber slides, or micro-dishes
3	<b>Imagery Wavelength</b>	Dyes: e.g. DAPI, FITC, Cy3, Cy5 Fluorescent proteins: e.g. EBFP2, EGFP, DsRed/mCherry
5	<b>Camera</b>	sCMOS camera (resolution: 2048 x 2048, pixel size: 6.5 μm x 6.5 μm)
	<b>Binning Options</b>	Low resolution mode: 800 x 800 pixels High resolution mode: 1600 x 1600 pixels
6	<b>Operating System</b>	Microsoft Windows 10
	<b>Pattern Segmentation Options</b>	Toolbox for traditional image processing Trained model using AI deep learning
<b>Labeling Resolution</b>		300 nm+**
<b>Sample Format</b>		Cells - fixed on a chambered coverslip Tissues - slide mounted FFPE (5 - 10 μm in thickness) or frozen tissue section (10 - 20 μm in thickness)
<b>Sample Size Requirement</b>		Cell numbers: 4 x 10 <sup>5</sup> to 1 x 10 <sup>6</sup> cells for a single LC-MS/MS analysis* Tissue slides: 4 - 8 tissue section for a single LC-MS/MS analysis*

\*Application dependent on and varying with the area and the number of ROIs. \*\*Objective dependent

All product specifications and data are subject to change without notice to improve reliability, function, design, or otherwise. For Research Use Only. Not for use in diagnostic procedures.



syncell.com / +1 617 631 2746  
 200 Dexter Ave, Watertown, MA 02472, USA  
 ©2024 Syncell, Inc. All rights reserved. Syncell, Microscop<sup>®</sup>, and Autoscoop<sup>™</sup> are pending trademarks in the United States and/or other countries. Doc#MC AN-003, Version 1.2.

